

## Product datasheet for **HP209452**

### WNT6 Human qPCR Primer Pair (NM\_006522)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7475
Forward Sequence:	GTGCAACTGCACAACAACGAGG
Reverse Sequence:	GAAATGGAGGCAGCTTCTGCCA
ACCN:	<u><a href="#">NM_006522</a></u> , <u><a href="#">NM_006522.1</a></u> , <u><a href="#">NM_006522.2</a></u> , <u><a href="#">NM_006522.3</a></u> , <u><a href="#">BC004329</a></u> , <u><a href="#">BC004329.1</a></u> , <u><a href="#">BT007456</a></u>
UniProt ID:	<u><a href="#">Q9Y6F9</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)