

Product datasheet for **HP209460**

ELL Human qPCR Primer Pair (NM_006532)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8178
Forward Sequence:	GTGCCATTGTCATCAAGGCTGG
Reverse Sequence:	CACCACTCTTCTGATGGCACT
ACCN:	BC010010 , NM_006532 , NM_006532.1 , NM_006532.2 , NM_006532.3 , BC049195 , BC049195.1 , BC033673 , BC064558 , BC108256 , BI458899 , NM_006532.4
UniProt ID:	P55199
Synonyms:	C19orf17; ELL1; MEN; PPP1R68
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)