

Product datasheet for **HP209474**

CAMKK2 Human qPCR Primer Pair (NM_006549)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10645
Forward Sequence:	CCTCAAACCACTCTCTGAAGACC
Reverse Sequence:	TTCTCCGACCAGGAGGTTGGAA
ACCN:	<u>NM_006549</u> , <u>NM_006549.1</u> , <u>NM_006549.2</u> , <u>NM_006549.3</u> , <u>BC026060</u> , <u>BC026060.2</u> , <u>BC000318</u> , <u>BF972593</u>
UniProt ID:	<u>Q96RR4</u>
Synonyms:	CAMKK; CAMKKB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)