

## Product datasheet for **HP209480**

### PMVK Human qPCR Primer Pair (NM\_006556)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10654
Forward Sequence:	GCCTTTCGGAAGGACATGATCC
Reverse Sequence:	ACTCTCCGTGTGCTACTACCA
ACCN:	<a href="#">NM_006556</a> , <a href="#">NM_006556.1</a> , <a href="#">NM_006556.2</a> , <a href="#">NM_006556.3</a> , <a href="#">BC006089</a> , <a href="#">BC006089.1</a> , <a href="#">BC007694</a> , <a href="#">BF221848</a> , <a href="#">BQ061044</a> , <a href="#">BT019976</a> , <a href="#">NM_006556.4</a>
UniProt ID:	<a href="#">Q15126</a>
Synonyms:	HUMPMKI; PMK; PMKA; PMKASE; POROK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)