

Product datasheet for **HP209492**

RRAGA Human qPCR Primer Pair (NM_006570)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10670
Forward Sequence:	TGAGACGCTCTACAAAGCCTGG
Reverse Sequence:	CGAACAGCAGAACTTCATCGGC
ACCN:	<u>NM_006570</u> , <u>NM_006570.3</u> , <u>NM_006570.4</u> , <u>BC009990</u> , <u>BC009990.2</u> , <u>BC006433</u> , <u>BX346838</u> , <u>BX442654</u> , <u>NM_006570.5</u>
UniProt ID:	<u>Q7L523</u>
Synonyms:	FIP-1; FIP1; RAGA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)