

Product datasheet for **HP209511**

C1orf2 (FAM189B) Human qPCR Primer Pair (NM_006589)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10712
Forward Sequence:	GTGGTTGGCATTGTGCCTGGA
Reverse Sequence:	TCCTTCCAGAGAGCACTGTTGG
ACCN:	<u>NM_006589</u> , <u>NM_006589.1</u> , <u>NM_006589.2</u> , <u>BC010471</u> , <u>BC010471.1</u> , <u>BC040136</u> , <u>BC040136.1</u> , <u>BC006297</u> , <u>BC006493</u> , <u>BC007444</u> , <u>BC008854</u> , <u>BC017262</u> , <u>BC033707</u> , <u>NM_006589.3</u>
UniProt ID:	<u>P81408</u>
Synonyms:	C1orf2; COTE1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)