

Product datasheet for **HP209513**

POLD3 Human qPCR Primer Pair (NM_006591)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10714
Forward Sequence:	TGAGCCTGAACCTCCTTCTGTC
Reverse Sequence:	CTGTGCAGGATTCCTCTCGTAG
ACCN:	<u>NM_006591</u> , <u>NM_006591.1</u> , <u>NM_006591.2</u> , <u>BC020587</u> , <u>BC032636</u> , <u>BC041703</u> , <u>BC058840</u> , <u>BC092454</u> , <u>BC108908</u> , <u>BC108909</u> , <u>NM_006591.3</u>
UniProt ID:	<u>Q15054</u>
Synonyms:	P66; P68; PPP1R128
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)