

Product datasheet for **HP209519**

KCC4 (SLC12A7) Human qPCR Primer Pair (NM_006598)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10723
Forward Sequence:	CCTCAAGGATGCACAGAAGTCC
Reverse Sequence:	CGTAAGACCACGCCTTCAATGC
ACCN:	BC007760 , NM_006598 , NM_006598.1 , NM_006598.2 , BC098390 , BC098390.1 , BC018982 , BU953746 , NM_006598.3
UniProt ID:	Q9Y666
Synonyms:	KCC4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)