

Product datasheet for **HP209527**

MASP2 Human qPCR Primer Pair (NM_006610)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10747
Forward Sequence:	GAGACCACACAGGCTGGAAGAT
Reverse Sequence:	GTCATAGCCAGTCTCGAAAAG
ACCN:	<u>NM_006610</u> , <u>NM_006610.1</u> , <u>NM_006610.2</u> , <u>NM_006610.3</u> , <u>BC052299</u> , <u>BC067359</u> , <u>BC080556</u> , <u>BC156086</u> , <u>BC156886</u> , <u>NM_006610.4</u>
UniProt ID:	<u>O00187</u>
Synonyms:	MAP-2; MAP19; MASP-2; MASP1P1; sMAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)