

## Product datasheet for **HP209543**

### ARPP19 Human qPCR Primer Pair (NM\_006628)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10776
Forward Sequence:	CAAAGCCTGGAGGTTTCAGATTTTC
Reverse Sequence:	GTCACCACTGACCTCCGTCTTA
ACCN:	<a href="#">NM_006628</a> , <a href="#">NM_006628.1</a> , <a href="#">NM_006628.2</a> , <a href="#">NM_006628.3</a> , <a href="#">NM_006628.4</a> , <a href="#">NM_006628.5</a> , <a href="#">BC100015</a> , <a href="#">BC100015.1</a> , <a href="#">BC003418</a> , <a href="#">BQ422616</a> , <a href="#">NM_006628.6</a>
UniProt ID:	<a href="#">P56211</a>
Synonyms:	ARPP-16; ARPP-19; ARPP16; ENSAL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)