

Product datasheet for **HP209565**

TUBGCP2 Human qPCR Primer Pair (NM_006659)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10844
Forward Sequence:	TACTCCGCTGTGACCAGGTTC A
Reverse Sequence:	CCAGAATCAGGTGCTCCTTCAC
ACCN:	<u>BC010373</u> , <u>NM_006659</u> , <u>NM_006659.1</u> , <u>NM_006659.2</u> , <u>NM_006659.3</u> , <u>BC093770</u> , <u>BC093770.1</u> , <u>BC005011</u> , <u>BC052368</u> , <u>BC111957</u> , <u>BC143247</u>
UniProt ID:	<u>Q9BSJ2</u>
Synonyms:	ALP4; GCP-2; GCP2; Grip103; h103p; hGCP2; hSpc97; PAMDDFS; SPBC97; SPC97; Spc97p
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)