

Product datasheet for **HP209566**

CLPX Human qPCR Primer Pair (NM_006660)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10845
Forward Sequence:	CAGGCTGGATATGTAGGCGAAG
Reverse Sequence:	TGAATGCCTGGCACACTGCCAA
ACCN:	<u>NM_006660</u> , <u>NM_006660.1</u> , <u>NM_006660.2</u> , <u>NM_006660.3</u> , <u>NM_006660.4</u> , <u>BC130373</u> , <u>BC136487</u> , <u>BC144166</u> , <u>BG717994</u> , <u>NM_006660.5</u>
UniProt ID:	<u>O76031</u>
Synonyms:	EPP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)