

Product datasheet for **HP209580**

USP19 Human qPCR Primer Pair (NM_006677)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10869
Forward Sequence:	GCTGCTATCCTCAGAGTTGGCT
Reverse Sequence:	TCATCCTCCGACTGTTGCTTCC
ACCN:	NM_006677 , NM_006677.1 , NM_006677.2 , BC014934 , BC048269 , BC065909 , BC082241 , BC106029 , BC128088 , BC142660 , BC142727 , BC146752
UniProt ID:	O94966
Synonyms:	ZMYND9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)