

Product datasheet for **HP209593**

JTB Human qPCR Primer Pair (NM_006694)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10899
Forward Sequence:	GGTGAAGAGTTTGTGGTAGCAG
Reverse Sequence:	CTCTTAGATGAGCTGCATGTGATT
ACCN:	NM_006694 , NM_006694.1 , NM_006694.2 , NM_006694.3 , BC000996 , BC000996.2 , BC000499 , BC001363 , BC001667 , BC004239 , BC019277 , BM907938 , BT007285 , BU527541 , NM_006694.4
UniProt ID:	O76095
Synonyms:	hJT; HJTB; HSPC222; PAR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)