

## Product datasheet for **HP209610**

### FASTK Human qPCR Primer Pair (NM\_006712)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10922
Forward Sequence:	AGGCAGCAAGTGCCCATCTTTC
Reverse Sequence:	ACAGTCAGGTCCTGGCGGTATT
ACCN:	<a href="#">NM_006712</a> , <a href="#">NM_006712.1</a> , <a href="#">NM_006712.2</a> , <a href="#">NM_006712.3</a> , <a href="#">NM_006712.4</a> , <a href="#">BC011770</a> , <a href="#">BC011770.2</a> , <a href="#">BC039026</a> , <a href="#">BC039026.2</a> , <a href="#">BC000377</a> , <a href="#">BC006386</a> , <a href="#">BM799420</a> , <a href="#">BU682096</a> , <a href="#">NM_006712.5</a>
UniProt ID:	<a href="#">Q14296</a>
Synonyms:	FAST
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)