

Product datasheet for **HP209613**

MAN2C1 Human qPCR Primer Pair (NM_006715)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4123
Forward Sequence:	CTCCTCATCGTCAACACACTGC
Reverse Sequence:	GAGGAACAGGAGCATAGCCCAT
ACCN:	<u>BC010081</u> , <u>NM_006715</u> , <u>NM_006715.1</u> , <u>NM_006715.2</u> , <u>NM_006715.3</u> , <u>BC050550</u> , <u>BC050550.1</u> , <u>BC038594</u> , <u>BC080191</u> , <u>BM151712</u> , <u>BX647111</u> , <u>NM_006715.4</u>
UniProt ID:	<u>Q9NTJ4</u>
Synonyms:	MAN6A8; MANA; MANA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)