

## Product datasheet for **HP209621**

### HIVEP2 Human qPCR Primer Pair (NM\_006734)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3097
Forward Sequence:	CGGTGGATGATACAGAACTGAG
Reverse Sequence:	CTGATTCCTCAGCATCGGAGAAC
ACCN:	<u><a href="#">NM_006734</a></u> , <u><a href="#">NM_006734.1</a></u> , <u><a href="#">NM_006734.2</a></u> , <u><a href="#">NM_006734.3</a></u> , <u><a href="#">BC032679</a></u> , <u><a href="#">BC167801</a></u> , <u><a href="#">NM_006734.4</a></u>
UniProt ID:	<u><a href="#">P31629</a></u>
Synonyms:	HIV-EP2; MBP-2; MIBP1; MRD43; SHN2; ZAS2; ZNF40B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)