

Product datasheet for **HP209626**

Protein phosphatase inhibitor 1 (PPP1R1A) Human qPCR Primer Pair (NM_006741)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5502
Forward Sequence:	CACAGAAGTGGAGTCAAGGCTG
Reverse Sequence:	TTGGCTCCCTTGGAAATCCAGTG
ACCN:	<u>NM_006741</u> , <u>NM_006741.1</u> , <u>NM_006741.2</u> , <u>NM_006741.3</u> , <u>BC022470</u> , <u>BC022470.1</u> , <u>BM876308</u> , <u>BP211831</u>
UniProt ID:	<u>Q13522</u>
Synonyms:	I1; IPP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)