

Product datasheet for **HP209631**

Syntrophin (SNTB2) Human qPCR Primer Pair (NM_006750)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6645
Forward Sequence:	CCACAGGTGTTGGCTGAACTCA
Reverse Sequence:	TGAGGACAGGTCTCCATTGCTG
ACCN:	<u>NM_006750</u> , <u>NM_006750.1</u> , <u>NM_006750.2</u> , <u>NM_006750.3</u> , <u>BC048215</u> , <u>BC031058</u> , <u>BC035134</u> , <u>BC036429</u>
UniProt ID:	<u>Q13425</u>
Synonyms:	D16S2531E; EST25263; SNT2B2; SNT3; SNTL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)