

## Product datasheet for **HP209683**

### **TOMM34 Human qPCR Primer Pair (NM\_006809)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	10953
<b>Forward Sequence:</b>	GTAGTAACCTGGAATCTGCCACG
<b>Reverse Sequence:</b>	CCGTCTGTAGAATGCCTTCACG
<b>ACCN:</b>	<a href="#">NM_006809</a> , <a href="#">NM_006809.1</a> , <a href="#">NM_006809.2</a> , <a href="#">NM_006809.3</a> , <a href="#">NM_006809.4</a> , <a href="#">BC007423</a> , <a href="#">BC007423.2</a> , <a href="#">BC001763</a> , <a href="#">BC006573</a> , <a href="#">BC014907</a> , <a href="#">BI560064</a> , <a href="#">BT020008</a> , <a href="#">BT020009</a> , <a href="#">NM_006809.5</a>
<b>UniProt ID:</b>	<a href="#">Q15785</a>
<b>Synonyms:</b>	HTOM34P; TOM34; URCC3
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)