

Product datasheet for **HP209707**

Cyclin (CCNI) Human qPCR Primer Pair (NM_006835)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10983
Forward Sequence:	CTTCACACAGCCACACCATTGG
Reverse Sequence:	GTAGTTGCTTGTAAGGACTGCC
ACCN:	<u>NM_006835</u> , <u>NM_006835.1</u> , <u>NM_006835.2</u> , <u>BC004975</u> , <u>BC004975.1</u> , <u>BC000420</u> , <u>BT020019</u> , <u>NM_006835.3</u>
UniProt ID:	<u>Q14094</u>
Synonyms:	CCNI1; CYC1; CYI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)