

Product datasheet for **HP209709**

JAB1 (COPS5) Human qPCR Primer Pair (NM_006837)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10987
Forward Sequence:	CCAGGAACCATTTGTAGCAGTGG
Reverse Sequence:	GTCTGGTACTCAGAAGGTCCTTC
ACCN:	<u>NM_006837</u> , <u>NM_006837.1</u> , <u>NM_006837.2</u> , <u>BC001187</u> , <u>BC001187.1</u> , <u>BC001859</u> , <u>BC007272</u> , <u>BI752436</u> , <u>BX648542</u> , <u>NM_006837.3</u>
UniProt ID:	<u>Q92905</u>
Synonyms:	CSN5; JAB1; MOV-34; SGN5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)