

Product datasheet for **HP209713**

SDS Human qPCR Primer Pair (NM_006843)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10993
Forward Sequence:	ACCAGCGTCTACCTCAAGATGG
Reverse Sequence:	GCCGAGGAGCAGACAAAATGTG
ACCN:	<u>NM_006843</u> , <u>NM_006843.1</u> , <u>NM_006843.2</u> , <u>BC020750</u> , <u>BG534053</u> , <u>BG566257</u> , <u>BU740747</u> , <u>NM_006843.3</u>
UniProt ID:	<u>P20132</u>
Synonyms:	SDH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)