

## Product datasheet for **HP209732**

### **RAB31 Human qPCR Primer Pair (NM\_006868)**

#### **Product data:**

Product Type:	qPCR Primer Pairs
Gene ID:	11031
Forward Sequence:	GGAGCTGAAAGAACATGGTCCAG
Reverse Sequence:	GCACTTGTCTCAACCACGATGG
ACCN:	<u><a href="#">NM_006868</a></u> , <u><a href="#">NM_006868.1</a></u> , <u><a href="#">NM_006868.2</a></u> , <u><a href="#">NM_006868.3</a></u> , <u><a href="#">BC001148</a></u> , <u><a href="#">BC001148.1</a></u> , <u><a href="#">BT020027</a></u> , <u><a href="#">BT020028</a></u> , <u><a href="#">BU608212</a></u> , <u><a href="#">BU944638</a></u> , <u><a href="#">NM_006868.4</a></u>
UniProt ID:	<u><a href="#">Q13636</a></u>
Synonyms:	Rab22B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)