

Product datasheet for **HP209742**

TIS11D (ZFP36L2) Human qPCR Primer Pair (NM_006887)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	678
Forward Sequence:	GTACGGCGAAAAGTGCCAGTTC
Reverse Sequence:	AGCCGATGGTATGAAAGGTGCG
ACCN:	<u>NM_006887</u> , <u>NM_006887.1</u> , <u>NM_006887.2</u> , <u>NM_006887.3</u> , <u>NM_006887.4</u> , <u>BC005010</u> , <u>BC020369</u> , <u>BE222042</u> , <u>BI489986</u> , <u>BI548620</u> , <u>NM_006887.5</u>
UniProt ID:	<u>P47974</u>
Synonyms:	BRF2; ERF-2; ERF2; RNF162C; TIS11D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)