

Product datasheet for **HP209751**

HOXD3 Human qPCR Primer Pair (NM_006898)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3232
Forward Sequence:	GTCTCGACAGAACTCCAAGCAG
Reverse Sequence:	CGCTCGTGTATGCCGTGCGTA
ACCN:	<u>NM_006898</u> , <u>NM_006898.1</u> , <u>NM_006898.2</u> , <u>NM_006898.3</u> , <u>NM_006898.4</u> , <u>BC005124</u> , <u>BC005124.2</u> , <u>BC008789</u> , <u>NM_006898.5</u>
UniProt ID:	<u>P31249</u>
Synonyms:	Hox-4.1; HOX1D; HOX4; HOX4A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)