

Product datasheet for **HP209759**

RASGRF2 Human qPCR Primer Pair (NM_006909)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5924
Forward Sequence:	GATTGAGGAGCCAGATGCAAGC
Reverse Sequence:	CAAGACAACAGTGAAGGCGGCA
ACCN:	NM_006909 , NM_006909.1 , NM_006909.2 , BC041953 , BC126112 , BC136296 , BC143870 , BU430075 , BX502901 , NM_006909.3
UniProt ID:	O14827
Synonyms:	GRF2; RAS-GRF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)