

Product datasheet for **HP209788**

SRP72 Human qPCR Primer Pair (NM_006947)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6731
Forward Sequence:	ATGGCTGCCAATGCGAGAACGT
Reverse Sequence:	GGCTGCTCACAGTTTACTGGC
ACCN:	<u>NM_006947</u> , <u>NM_006947.1</u> , <u>NM_006947.2</u> , <u>NM_006947.3</u> , <u>BC017057</u> , <u>BC032609</u> , <u>BC040134</u> , <u>BC046143</u> , <u>BC056901</u> , <u>BC065018</u> , <u>BC105583</u> , <u>BC160164</u> , <u>BM985263</u> , <u>BM985419</u> , <u>BX537991</u> , <u>NM_006947.4</u>
UniProt ID:	<u>O76094</u>
Synonyms:	BMFF; BMFS1; HEL103
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)