

## Product datasheet for **HP209794**

### Uroplakin III (UPK3A) Human qPCR Primer Pair (NM\_006953)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7380
Forward Sequence:	CTCACAGATCCTGAATGCCTACC
Reverse Sequence:	CCGTGGACATATTGACCAGGAC
ACCN:	<a href="#">NM_006953</a> , <a href="#">NM_006953.1</a> , <a href="#">NM_006953.2</a> , <a href="#">NM_006953.3</a> , <a href="#">BC069544</a> , <a href="#">BC069544.1</a> , <a href="#">BC108900</a> , <a href="#">BM681088</a> , <a href="#">BM727658</a> , <a href="#">NM_006953.4</a>
UniProt ID:	<a href="#">O75631</a>
Synonyms:	UP3A; UPIII; UPIIIA; UPK3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)