

## Product datasheet for **HP209821**

### TRF41 (PAPD7) Human qPCR Primer Pair (NM\_006999)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11044
Forward Sequence:	CAATGCCAGTGGCAAACCTCA
Reverse Sequence:	ACGCAGTTCCTTCTACACCAGC
ACCN:	<u><a href="#">NM_006999</a></u> , <u><a href="#">NM_006999.1</a></u> , <u><a href="#">NM_006999.2</a></u> , <u><a href="#">NM_006999.3</a></u> , <u><a href="#">NM_006999.4</a></u> , <u><a href="#">BC084567</a></u> , <u><a href="#">BC117137</a></u> , <u><a href="#">BC126106</a></u> , <u><a href="#">BC143882</a></u>
UniProt ID:	<u><a href="#">Q5XG87</a></u>
Synonyms:	LAK-1; LAK1; PAPD7; POLK; POLS; TRF4; TRF4-1; TRF41; TUTASE5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)