

Product datasheet for **HP209825**

PAGE4 Human qPCR Primer Pair (NM_007003)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9506
Forward Sequence:	CGGTGAATCTCAGCAAGAGGAAC
Reverse Sequence:	GCTCACTCCGAGTCTTTTCCAG
ACCN:	<u>NM_007003</u> , <u>NM_007003.2</u> , <u>NM_007003.3</u> , <u>NM_007003.4</u> , <u>BC010897</u> , <u>BC010897.2</u> , <u>BG435968</u> , <u>BG622311</u> , <u>BX371732</u>
UniProt ID:	<u>O60829</u>
Synonyms:	CT16.7; GAGE-9; GAGEC1; JM-27; JM27; PAGE-1; PAGE-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)