

## Product datasheet for **HP209832**

### ABHD2 Human qPCR Primer Pair (NM\_007011)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11057
Forward Sequence:	CTGACCTCTACTTCCAGGACTC
Reverse Sequence:	CCATACAAGGCTGTCTGGATGTG
ACCN:	<u><a href="#">NM_007011</a></u> , <u><a href="#">NM_007011.1</a></u> , <u><a href="#">NM_007011.2</a></u> , <u><a href="#">NM_007011.3</a></u> , <u><a href="#">NM_007011.5</a></u> , <u><a href="#">NM_007011.6</a></u> , <u><a href="#">NM_007011.7</a></u> , <u><a href="#">BC019248</a></u> , <u><a href="#">BC019248.1</a></u> , <u><a href="#">BC090052</a></u> , <u><a href="#">BC090052.1</a></u> , <u><a href="#">BC017905</a></u> , <u><a href="#">BI459790</a></u> , <u><a href="#">BI464471</a></u>
UniProt ID:	<u><a href="#">P08910</a></u>
Synonyms:	HS1-2; LABH2; PPHP1-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)