

## Product datasheet for **HP209853**

### ADAMTS8 Human qPCR Primer Pair (NM\_007037)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11095
Forward Sequence:	GCTGTCTACCTGAGGAGGAAGT
Reverse Sequence:	TCCTCCACAGGTCCGAGAACAT
ACCN:	<a href="#">NM_007037</a> , <a href="#">NM_007037.1</a> , <a href="#">NM_007037.2</a> , <a href="#">NM_007037.3</a> , <a href="#">NM_007037.4</a> , <a href="#">NM_007037.5</a> , <a href="#">BC089435</a> , <a href="#">BC113874</a> , <a href="#">BC114624</a> , <a href="#">BM988619</a> , <a href="#">BX424609</a> , <a href="#">NM_007037.6</a>
UniProt ID:	<a href="#">Q9UP79</a>
Synonyms:	ADAM-TS8; METH2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)