

Product datasheet for **HP209854**

ADAMTS5 Human qPCR Primer Pair (NM_007038)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11096
Forward Sequence:	CCTGGTCCAAATGCACTTCAGC
Reverse Sequence:	TCGTAGGTCTGTCCTGGGAGTT
ACCN:	NM_007038 , NM_007038.1 , NM_007038.2 , NM_007038.3 , NM_007038.4 , BC093775 , BC093775.1 , BC093777
UniProt ID:	Q9UNA0
Synonyms:	ADAM-TS 5; ADAM-TS5; ADAM-TS 11; ADAMTS-5; ADAMTS-11; ADAMTS11; ADMP-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)