

Product datasheet for **HP209874**

PWP1 Human qPCR Primer Pair (NM_007062)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11137
Forward Sequence:	TAGAGCCAGTCTTCACACTCGG
Reverse Sequence:	GGTCAAGGACAGCATCGGTATG
ACCN:	NM_007062 , NM_007062.1 , NM_007062.2 , BC001652 , BC001652.2 , BC000067 , BC010921 , BC032127 , BC040135 , BC046144 , NM_007062.3
UniProt ID:	Q13610
Synonyms:	IEF-SSP-9502
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)