

Product datasheet for **HP209878**

DMC1 Human qPCR Primer Pair (NM_007068)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11144
Forward Sequence:	AATGACTGCCGATCCAGGAGCA
Reverse Sequence:	TCTCCTCTCCCTTCGCAAGC
ACCN:	<u>NM_007068</u> , <u>NM_007068.1</u> , <u>NM_007068.2</u> , <u>NM_007068.3</u> , <u>BC125164</u> , <u>BC035658</u> , <u>BC125163</u> , <u>BG959880</u> , <u>BM545092</u> , <u>BM843634</u> , <u>NM_007068.4</u>
UniProt ID:	<u>Q14565</u>
Synonyms:	dj199H16.1; DMC1H; LIM15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)