

Product datasheet for **HP209883**

Coronin 1a (CORO1A) Human qPCR Primer Pair (NM_007074)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11151
Forward Sequence:	CCAACATCGTCTACCTCTGTGG
Reverse Sequence:	CTCACACTTGTTACCTCCAGG
ACCN:	NM_007074 , NM_007074.1 , NM_007074.2 , NM_007074.3 , BC110374 , BC110374.1 , BC126385 , BC126387
UniProt ID:	P31146
Synonyms:	CLABP; CLIPINA; HCORO1; IMD8; p57; TACO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)