

Product datasheet for **HP209892**

ATP5ME Human qPCR Primer Pair (NM_007100)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	521
Forward Sequence:	CCTAAAACCTCGGGCAGAAGAG
Reverse Sequence:	GGTCACTCACTTTAATATGCTGTC
ACCN:	<u>NM_007100</u> , <u>NM_007100.1</u> , <u>NM_007100.2</u> , <u>NM_007100.3</u> , <u>BC105610</u> , <u>BC003679</u> , <u>BC017215</u> , <u>BF215596</u> , <u>BG940181</u> , <u>BG943252</u> , <u>BT006922</u> , <u>BT008004</u> , <u>NM_007100.4</u>
UniProt ID:	<u>P56385</u>
Synonyms:	ATP5I; ATP5K
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)