

Product datasheet for **HP209895**

NDUFV1 Human qPCR Primer Pair (NM_007103)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4723
Forward Sequence:	TGTGTGAGACGGTGCTGATGGA
Reverse Sequence:	CGATGGCTTTCACGATGTCCGT
ACCN:	<u>NM_007103</u> , <u>NM_007103.1</u> , <u>NM_007103.2</u> , <u>NM_007103.3</u> , <u>BC015645</u> , <u>BC015645.2</u> , <u>BC008146</u> , <u>BC007619</u> , <u>BP254989</u> , <u>NM_007103.4</u>
UniProt ID:	<u>P49821</u>
Synonyms:	CI-51K; CI51KD; MC1DN4; UQOR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)