

Product datasheet for **HP209897**

SLC22A18AS Human qPCR Primer Pair (NM_007105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5003
Forward Sequence:	TCCAGGCTGACTGGAGGAAGG
Reverse Sequence:	GAGAACCACCAAGCACCCCTG
ACCN:	<u>BC030237</u> , <u>NM_007105</u> , <u>NM_007105.1</u> , <u>NM_007105.2</u> , <u>NM_007105.3</u> , <u>BC030237.1</u> , <u>NM_007105.4</u>
UniProt ID:	<u>Q8N1D0</u>
Synonyms:	BWR1B; BWSCR1B; ORCTL2S; p27-BWR1B; SLC22A1LS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)