

Product datasheet for **HP209909**

UGT1A4 Human qPCR Primer Pair (NM_007120)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54657
Forward Sequence:	GCAAAGCGCATGGAGACTAAGG
Reverse Sequence:	GGTCCTTGTGAAGGCTGGAGAG
ACCN:	<u>NM_007120</u> , <u>NM_007120.1</u> , <u>NM_007120.2</u> , <u>BC121036</u> , <u>BC131623</u> , <u>BC139784</u> , <u>NM_007120.3</u>
UniProt ID:	<u>P22310</u>
Synonyms:	GNT1; hUG-BR1; HUG-BR2; UDPGT; UDPGT 1-4; UGT-1A; UGT-1D; UGT1; UGT1-01; UGT1-04; UGT1.1; UGT1.4; UGT1A; UGT1A1; UGT1A4S; UGT1D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)