

Product datasheet for **HP209911**

USF1 Human qPCR Primer Pair (NM_007122)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7391
Forward Sequence:	GCTCTATGGAGAGCACCAAGTC
Reverse Sequence:	AGACAAGCGGTGGTTACTCTGC
ACCN:	NM_007122 , NM_007122.1 , NM_007122.2 , NM_007122.3 , NM_007122.4 , BC035505 , BC035505.2 , BP280296 , NM_007122.5
UniProt ID:	P22415
Synonyms:	bHLHb11; FCHL; FCHL1; HYPLIP1; MLTF; MLTFI; UEF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)