

## Product datasheet for **HP209949**

### Unrip (STRAP) Human qPCR Primer Pair (NM\_007178)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11171
Forward Sequence:	GATGCCACCAAAGCAGCTACAG
Reverse Sequence:	CCTGCGTGAAATCCACAGTCTTG
ACCN:	<u><a href="#">NM_007178</a></u> , <u><a href="#">NM_007178.1</a></u> , <u><a href="#">NM_007178.2</a></u> , <u><a href="#">NM_007178.3</a></u> , <u><a href="#">BC062306</a></u> , <u><a href="#">BC062306.1</a></u> , <u><a href="#">BC000162</a></u> , <u><a href="#">BT020044</a></u> , <u><a href="#">BT020045</a></u> , <u><a href="#">NM_007178.4</a></u>
UniProt ID:	<u><a href="#">Q9Y3F4</a></u>
Synonyms:	MAWD; PT-WD; UNRIP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)