

Product datasheet for **HP209976**

RNF2 Human qPCR Primer Pair (NM_007212)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6045
Forward Sequence:	CAGTCACAGCATTGAGGAAGGAC
Reverse Sequence:	GCTTCCTGATTGCTATGTGTGGA
ACCN:	<u>NM_007212</u> , <u>NM_007212.1</u> , <u>NM_007212.2</u> , <u>NM_007212.3</u> , <u>BC012583</u> , <u>BC012583.1</u> , <u>BF183837</u> , <u>BX338357</u> , <u>NM_007212.4</u>
UniProt ID:	<u>Q99496</u>
Synonyms:	BAP-1; BAP1; DING; HIPI3; RING1B; RING2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)