

Product datasheet for **HP209981**

PDCD10 Human qPCR Primer Pair (NM_007217)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11235
Forward Sequence:	CAGACACTGAGAGCCGCTTTCA
Reverse Sequence:	CTGCCATACGAAGAAGGGACTC
Accession No:	<u>NM_007217</u> , <u>NM_145859</u> , <u>NM_145860</u> , <u>BC002506</u> , <u>BC016353</u> , <u>NM_007217.4</u>
Synonyms:	CCM3; TFAR15
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)