

## Product datasheet for **HP209982**

### RNF139 Human qPCR Primer Pair (NM\_007218)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11236
Forward Sequence:	TGCATCGTGCTCCAGATCTTCC
Reverse Sequence:	CAGAAAGGCTGAGCTGTACGTG
ACCN:	<a href="#">NM_007218</a> , <a href="#">NM_007218.1</a> , <a href="#">NM_007218.2</a> , <a href="#">NM_007218.3</a> , <a href="#">BC064636</a> , <a href="#">BC021571</a> , <a href="#">BI562637</a> , <a href="#">BM677460</a> , <a href="#">BM804626</a> , <a href="#">NM_007218.4</a>
UniProt ID:	<a href="#">Q8WU17</a>
Synonyms:	HRCA1; RCA1; TRC8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)