

## Product datasheet for **HP210005**

### SYNRG Human qPCR Primer Pair (NM\_007247)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11276
Forward Sequence:	AGTGAGCTGGCAACCTTTGACC
Reverse Sequence:	TGGAACGGTCTCTGAAAGGCTG
ACCN:	<a href="#">NM_007247</a> , <a href="#">NM_007247.1</a> , <a href="#">NM_007247.2</a> , <a href="#">NM_007247.4</a> , <a href="#">NM_007247.5</a> , <a href="#">BC143476</a> , <a href="#">BC143478</a> , <a href="#">BC117313</a> , <a href="#">BC084536</a> , <a href="#">BC090930</a> , <a href="#">BC143452</a> , <a href="#">BG120023</a> , <a href="#">NM_007247.6</a>
UniProt ID:	<a href="#">Q9UMZ2</a>
Synonyms:	AP1GBP1; SYNG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)