

Product datasheet for **HP210023**

FKBP9 Human qPCR Primer Pair (NM_007270)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11328
Forward Sequence:	GAAAAGCGAAGGATTGTGGTCCC
Reverse Sequence:	TGCTGATGGAGTCCGAAGGGTT
ACCN:	<u>BD123495</u> , <u>NM_007270</u> , <u>NM_007270.1</u> , <u>NM_007270.2</u> , <u>NM_007270.3</u> , <u>NM_007270.4</u> , <u>BC007443</u> , <u>BC064418</u> , <u>BC072422</u> , <u>BC099918</u> , <u>BC101723</u> , <u>BC112053</u> , <u>BC127737</u> , <u>BC150541</u> , <u>BC150542</u> , <u>BQ688269</u> , <u>NM_007270.5</u>
UniProt ID:	<u>O95302</u>
Synonyms:	FKBP60; FKBP63; PPIase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)