

Product datasheet for **HP210029**

HP1 gamma (CBX3) Human qPCR Primer Pair (NM_007276)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11335
Forward Sequence:	GCTGACAAACCAAGAGGATTTGC
Reverse Sequence:	CAGCACCAAGTCTGCCTCATCT
ACCN:	NM_007276 , NM_007276.2 , NM_007276.3 , NM_007276.4 , BC000954 , BC000954.1 , BC014380 , BX647444 , BX648807 , NM_007276.5
UniProt ID:	Q13185
Synonyms:	HECH; HP1-GAMMA; HP1Hs-gamma
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)