

## Product datasheet for **HP210048**

### ATF6 Human qPCR Primer Pair (NM\_007348)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 22926   |
| Forward Sequence: | CAGACAGTACCAACGCTTATGCC   |
| Reverse Sequence: | GCAGAACTCCAGGTGCTTGAAG  |
| ACCN:             | <u><a href="#">NM_007348</a></u> , <u><a href="#">NM_007348.1</a></u> , <u><a href="#">NM_007348.2</a></u> , <u><a href="#">NM_007348.3</a></u> , <u><a href="#">BC014969</a></u> , <u><a href="#">BC071997</a></u> , <u><a href="#">BC166676</a></u> ,<br><u><a href="#">BM683233</a></u> , <u><a href="#">BP380917</a></u> , <u><a href="#">BX538263</a></u> , <u><a href="#">NM_007348.4</a></u> |
| UniProt ID:       | <u><a href="#">P18850</a></u>   |
| Synonyms:         | ACHM7; ATF6A  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)